



Dumont #5 Mirror Finish Forceps - Biology Tips/Straight/Inox/11cm



These forceps are highly polished to prevent embryos, nerves, or other adhesive tissues from sticking to the tips.

General Info

Product Number	11252-23
Category	Forceps & Tweezers > Smooth
EAN	4251909502225



Find 11252-23 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	11 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Style	#5		
Tipdimensions	0.05 x 0.02 mm		
Tips	Biology		
Tipshape	Straight		