



Dumont #5 Forceps - Biology Tips/Straight/Dumoxel/11cm



General Info

Product Number	11252-30
Category	Forceps & Tweezers > Smooth Most Popular > Forceps & Tweezers
EAN	4251909502232



Find 11252-30 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Dumoxel		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Length	11 cm		
Style	#5		
Tipdimensions	0.05 x 0.02 mm		
Tips	Biology		
Tipshape	Straight		