



General Info

Product Number	11252-40
Category	Forceps & Tweezers > Smooth MedTech Line Most Popular > Forceps & Tweezers
EAN	4251909502249



Find 11252-40 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	7 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Titanium		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Length	11 cm		
Style	#5		
Tipdimensions	0.05 x 0.02 mm		
Tips	Biology		
Tipshape	Straight		