



Dumont #5 Forceps - Biology Tips/Straight/Ceramic Coated Inox/11cm



General Info

Product Number	11252-50
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	4251909502256



Find 11252-50 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Ceramic coated Surgical Stainless Steel (Inox)		
Length	11 cm		
Style	#5		
Tipdimensions	0.05 x 0.02 mm		
Tips	Biology		
Tipshape	Straight		