



Dumont #5XL Forceps - Standard Tips/Straight/Inox/15cm



General Info

Product Number	11253-10
Category	Forceps & Tweezers > Smooth
EAN	7640241754671



Find 11253-10 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	20 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15 cm		
Style	#5XL		
Tipdimensions	0.1 x 0.06 mm		
Tips	Standard		
Tipshape	Straight		