



Dumont #5A Medical Forceps - Micro-Blunted Standard Tips/Off Center/Inox/11cm



These forceps have atraumatic or gently blunted tips to minimize the chance of puncturing tissue.

General Info

Product Number	11253-21
Category	Forceps & Tweezers > Smooth Most Popular > Forceps & Tweezers
EAN	7640241754657



Find 11253-21 on finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	14 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11.5 cm		
Style	Medical #5A		
Tipdimensions	0.13 x 0.08 mm		
Tips	Standard		
Tipshape	Off-Center		