



Dumont #L5 Forceps -Sliding Lock/Standard Tips/Straight/Inox/11cm



Plastic clamping ring holds tips together and can be locked in place. Fully autoclavable.

General Info

Product Number	11253-27
Category	Forceps & Tweezers > Smooth
EAN	7640241754725



Find 11253-27 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Style	#L5		
Tipdimensions	0.1 x 0.06 mm		
Tips	Standard		
Tipshape	Straight		