



Dumont #5/90 Medical Forceps - Micro-Blunted Standard Tips/Angled
90°/Inox/11cm



These forceps have atraumatic or gently blunted tips to minimize the chance of puncturing tissue.

General Info

Product Number	11253-29
Category	Forceps & Tweezers > Smooth Most Popular > Forceps & Tweezers
EAN	7640241754640



Find 11253-29 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Style	Medical #5/90		
Tipdimensions	0.1 x 0.06 mm		
Tips	Standard		
Tipshape	Angled 90°		