



Dumont #55 Forceps - Biology Tips/Straight/Inox/11cm



Beveled like the #5 but with lighter and finer shanks.



General Info

Product Number	11255-20
Category	Forceps & Tweezers > Smooth Most Popular > Forceps & Tweezers
EAN	4251909502270



Find 11255-20 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Style	#55		
Tipdimensions	0.05 x 0.02 mm		
Tips	Biology		
Tipshape	Straight		