



Dumont #6 Forceps - Standard Tips/Angled/Inox/11.5cm



Sharply angled tips for access to difficult areas.

General Info

Product Number	11260-20
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	4251909502294



Find 11260-20 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	14 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11.5 cm		
Style	#6		
Tipdimensions	0.17 x 0.1 mm		
Tips	Standard		
Tipshape	Angled		