



Dumont #7B Forceps - Standard Tips/Curved/Inox/11cm



Allgemein

Produktnummer	11270-20
Kategorie	Forceps & Tweezers > Serrated
EAN	7640241754695



Finde 11270-20 auf [finescience.com](https://www.finescience.com)

Technische Informationen

Versandgewicht	15 g	Nettogewicht	14 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11.5 cm		
Serrated	Yes		
Style	#7b		
Tipdimensions	0.17 x 0.1 mm		
Tips	Standard		
Tipshape	Curved		