



Dumont #7 Forceps - Standard Tips/Curved/Dumoxel/11.5cm



Allgemein

Produktnummer	11271-30
Kategorie	Forceps & Tweezers > Smooth Most Popular > Forceps & Tweezers
EAN	7640241754381



Finde 11271-30 auf finescience.com

Technische Informationen

Versandgewicht	15 g	Nettogewicht	14 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Dumoxel		
Antimagnetic	Yes		
Length	11.5 cm		
Style	#7		
Tipdimensions	0.17 x 0.1 mm		
Tips	Standard		
Tipshape	Curved		