



Dumont #7 Forceps - Standard Tips/Curved/Dumoxel/11.5cm



General Info

Product Number	11271-30
Category	Forceps & Tweezers > Smooth Most Popular > Forceps & Tweezers
EAN	7640241754381



Find 11271-30 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	14 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Dumoxel		
Antimagnetic	Yes		
Length	11.5 cm		
Style	#7		
Tipdimensions	0.17 x 0.1 mm		
Tips	Standard		
Tipshape	Curved		