



Dumont #7 Forceps - Biology Tips/Curved/Dumoxel/11.5cm



General Info

Product Number	11272-30
Category	Forceps & Tweezers > Smooth Most Popular > Forceps & Tweezers
EAN	4251909502300



Find 11272-30 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	14 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Dumoxel		
Antimagnetic	Yes		
Length	11.5 cm		
Style	#7		
Tipdimensions	0.07 x 0.04 mm		
Tips	Biology		
Tipshape	Curved		