



Dumont #7 Forceps - Biology Tips/Curved/Ceramic Coated Inox/11.5cm



General Info

Product Number	11272-50
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	4251909502324



Find 11272-50 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	15 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Ceramic coated Surgical Stainless Steel (Inox)		
Length	11.5 cm		
Style	#7		
Tipdimensions	0.07 x 0.04 mm		
Tips	Biology		
Tipshape	Curved		