



Dumont #7 Medical Forceps - Micro-Blunted Standard
Tips/Curved/Inox/11.5cm



These forceps have atraumatic or gently blunted tips to minimize the chance of puncturing tissue.

General Info

Product Number	11273-20
Category	Forceps & Tweezers > Smooth Most Popular > Forceps & Tweezers
EAN	7640241754688



Find 11273-20 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	15 g	Net Weight	14 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11.5 cm		
Style	Medical #7		
Tipdimensions	0.17 x 0.1 mm		
Tips	Standard		
Tipshape	Curved		