



Dumont #7 Medical Forceps - Biology Tips/Curved/Inox/11.5cm



Finished under magnification, the tips close flat to 3 mm and therefore less likely to open under excess pressure.

General Info

Product Number	11274-20
Category	Forceps & Tweezers > Smooth Most Popular > Forceps & Tweezers
EAN	4251909502331



Find 11274-20 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	14 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11.5 cm		
Style	#7		
Tipdimensions	0.07 x 0.03 mm		
Tips	Biology		
Tipshape	Curved		