



Dumont Vessel Cannulation Forceps - 0.5mm OD/13cm



These vessel cannulation forceps are useful for introducing a small plastic tube into a blood vessel of almost equal size. The hollowed jaws hold the tubing securely without deforming it. The tip of the tubing is guided into the vessel opening while preventing damage from unwanted movement. For tubing of 0.5mm OD.

General Info

Product Number	11281-11
Category	Forceps & Tweezers > Cannulation
EAN	4251909502355



Find 11281-11 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	20 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Tipshape	Angled		
Tubingsize	0.5 mm OD		