



Dumont #2A Mini Forceps - Standard Tips/Straight/Dumostar/7cm



General Info

Product Number	11290-21
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	4251909511944



Find 11290-21 on
finescience.com

Technical Info

Shipping Weight	9 g	Net Weight	6 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Dumostar		
Antimagnetic	Yes		
Length	7 cm		
Style	Mini #2A		
Tipdimensions	1.2 x 0.09 mm		
Tips	Standard		
Tipshape	Straight		