



Dumont #3 Forceps - Standard Tips/Straight/Dumostar/12cm



Allgemein

Produktnummer	11293-00
Kategorie	Forceps & Tweezers > Smooth
EAN	4251909502379



Finde 11293-00 auf finescience.com

Technische Informationen

Versandgewicht	15 g	Nettogewicht	11 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Dumostar		
Antimagnetic	Yes		
Length	12 cm		
Style	#3		
Tipdimensions	0.17 x 0.1 mm		
Tips	Standard		
Tipshape	Straight		