



Dumont #4 Forceps - Standard Tips/Straight/Dumostar/11cm



Similar to #5 but with broader shanks.

General Info

Product Number	11294-00
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	4251909502386



Find 11294-00 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Dumostar		
Antimagnetic	Yes		
Length	11 cm		
Style	#4		
Tipdimensions	0.13 x 0.08 mm		
Tips	Standard		
Tipshape	Straight		