

F · S · T

Dumont #5CO Forceps - Biology Tips/Straight/Dumostar/11cm



Allgemein

Produktnummer	11295-20
Kategorie	Forceps & Tweezers > Smooth
EAN	4251909502416



Finde 11295-20 auf
[finescience.com](https://www.finescience.com)

Technische Informationen

Versandgewicht	15 g	Nettogewicht	11 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Dumostar		
Autoclavesafe	Yes		
Length	11 cm		
Style	#5CO		
Tipdimensions	0.05 x 0.02 mm		
Tips	Biology		
Tipshape	Straight		