



Dumont #5CO Forceps - Biology Tips/Straight/Dumostar/11cm



General Info

Product Number	11295-20
Category	Forceps & Tweezers > Smooth
EAN	4251909502416



Find 11295-20 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	11 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Dumostar		
Antimagnetic	Yes		
Length	11 cm		
Style	#5CO		
Tipdimensions	0.05 x 0.02 mm		
Tips	Biology		
Tipshape	Straight		