



#2AF Medical Forceps - Standard Tips/Straight/Superresistant Alloy/12cm



Allgemein

|               |                                          |
|---------------|------------------------------------------|
| Produktnummer | 11320-21                                 |
| Kategorie     | Forceps & Tweezers > Smooth MedTech Line |
| EAN           | 4251909511883                            |



Finde 11320-21 auf [finescience.com](https://finescience.com)

Technische Informationen

|                     |                                                   |              |      |
|---------------------|---------------------------------------------------|--------------|------|
| Versandgewicht      | 20 g                                              | Nettogewicht | 14 g |
| Paketgröße          | 1 Pc.                                             | Maßeinheit   | EA   |
| Mindestbestellmenge | 1                                                 |              |      |
| Alloymaterial       | Anti-Magnetic Anti-Acid Superalloy (CX), polished |              |      |
| Antimagnetic        | Yes                                               |              |      |
| Length              | 12 cm                                             |              |      |
| Style               | #2AF                                              |              |      |
| Tipdimensions       | 1.5 x 0.2 mm                                      |              |      |
| Tips                | Blunt                                             |              |      |
| Tipshape            | Straight                                          |              |      |