

F · S · T

Moria MC31 Iris Forceps - Serrated/Curved/10cm



Balanced with gentle spring action and tips that are perfectly matched.

General Info

Product Number	11370-31
Category	Forceps & Tweezers > Serrated
EAN	4251909502454



Find 11370-31 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	30 g	Net Weight	7 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloy material	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	10 cm		
Tip dimensions	0.5 x 0.3 mm		
Tip width	0.5 mm		
Tips	Serrated		
Tip shape	Curved		