



Moria MC40B #5 Forceps - Biology Tips/Curved/11cm



These forceps have been expertly finished under 40x magnification and have the most finely executed points available. They are extremely light, superbly balanced, and have a gentle spring action. The shanks are specially designed with an arch which ensures maximum force to the tips. The tips close flat for 6mm from the points to give a large holding area, and the bottom portion is slim to reduce their weight.

General Info

Product Number	11370-42
Category	Forceps & Tweezers > Smooth
EAN	4251909502485



Find 11370-42 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	14 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Tipdimensions	0.05 x 0.01 mm		
Tips	Biology		
Tipshape	Curved		