



Moria MC32B Iris Forceps - Smooth/Straight/10cm



Balanced with gentle spring action and tips that are perfectly matched.

General Info

Product Number	11373-22
Category	Forceps & Tweezers > Smooth
EAN	4251909502508



Find 11373-22 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	7 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	10 cm		
Tipdimensions	0.5 x 0.3 mm		
Tipwidth	0.5 mm		
Tips	Smooth		
Tipshape	Straight		