

F · S · T

## Moria 7835 Bonn Forceps - Tying Platform with 1x2 Teeth/Straight/Stainless Steel/7.5cm



### Allgemein

|               |   |
|---------------|---|
| Produktnummer | 11383-50  |
| Kategorie     | Forceps & Tweezers > Toothed<br>Forceps & Tweezers > Smooth |
| EAN           | 4251909512323   |



Finde 11383-50 auf  
[finescience.com](https://www.finescience.com)

### Technische Informationen

|                     |                                 |              |     |
|---------------------|---------------------------------|--------------|-----|
| Versandgewicht      | 26 g                            | Nettogewicht | 6 g |
| Paketgröße          | 1 Pc.                           | Maßeinheit   | EA  |
| Mindestbestellmenge | 1                               |              |     |
| Alloymaterial       | Surgical Stainless Steel (Inox) |              |     |
| Autoclavesafe       | Yes                             |              |     |
| Length              | 7.5 cm                          |              |     |
| Tipdimensions       | 0.3 x 0.4 mm                    |              |     |
| Tips                | Tying platform with 1x2 teeth   |              |     |
| Tipshape            | Straight                        |              |     |
| Tyingplatform       | 4 mm                            |              |     |