



Moria TI7835 Bonn Forceps - Tying Platform with 1x2
Teeth/Straight/Titanium/7.5cm



General Info

Product Number	11383-56
Category	Forceps & Tweezers > Smooth Forceps & Tweezers > Toothed
EAN	4251909512330



Find 11383-56 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	24 g	Net Weight	4 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Titanium		
Length	7.5 cm		
Tipdimensions	0.3 x 0.4 mm		
Tips	Tying platform with 1x2 teeth		
Tipshape	Straight		
Tyingplatform	4 mm		