

F · S · T

## Moria 3225 Paufique Forceps - Tying Platform with 1x2 Teeth/Straight/Stainless Steel/8.5cm



### Allgemein

Produktnummer	11390-25
Kategorie	Forceps & Tweezers > Toothed Forceps & Tweezers > Smooth
EAN	4251909512347



Finde 11390-25 auf  
[finescience.com](https://www.finescience.com)

### Technische Informationen

Versandgewicht	28 g	Nettogewicht	8 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	8.5 cm		
Tipdimensions	0.6 x 1 mm		
Tips	Tying platform with 1x2 teeth		
Tipshape	Straight		
Tyingplatform	5 mm		