



Moria 9604 Troutman Suture Tying Forceps - 9.5 cm/Straight/8 mm



Used to suture after surgery is complete. This requires the manipulation of a calibrated nylon thread and the use of a specific instrument: a tying forceps, straight or curved, without teeth.

General Info

Product Number	11393-04
Category	Forceps & Tweezers > Smooth
EAN	4251909512248



Find 11393-04 on finescience.com

Technical Info

Shipping Weight	27 g	Net Weight	7 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	9.5 cm		
Tipdiameter	0.3 mm		
Tipdimensions	0.3 x 0.4 mm		
Tips	Smooth with platform		
Tipshape	Straight		
Tyingplatform	8 mm		