



Moria 9605 Troutman Suture Tying Forceps - 8.5 cm/Curved/7.5 mm



Used to suture after surgery is complete. This requires the manipulation of a calibrated nylon thread and the use of a specific instrument: a tying forceps, straight or curved, without teeth.

General Info

Product Number	11393-05
Category	Forceps & Tweezers > Smooth
EAN	4251909512255



Find 11393-05 on finescience.com

Technical Info

Shipping Weight	29 g	Net Weight	9 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	8.5 cm		
Tipdiameter	0.3 mm		
Tipdimensions	0.3 x 0.4 mm		
Tips	Smooth with platform		
Tipshape	Curved		
Tyingplatform	7.5 mm		