



Moria 9987 #5XL Forceps - Standard Tips/Curved/13cm



These forceps have been expertly finished under 40x magnification and have the most finely executed points available. They are extremely light, superbly balanced, and have a gentle spring action. The shanks are specially designed with an arch which ensures maximum force to the tips. The tips close flat for 6mm from the points to give a large holding area. The bottom portion is slim to reduce their weight, and the fenestrated handle provides a secure grip and reduces their weight as well.

General Info

Product Number	11399-87
Category	Forceps & Tweezers > Smooth
EAN	4251909502522



Find 11399-87 on
finescience.com

Technical Info

Shipping Weight	55 g	Net Weight	17 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	13 cm		
Tipdimensions	0.1 x 0.06 mm		
Tips	Standard		
Tipshape	Curved		