



#5 Mirror Finish Forceps - Standard Tips/Straight/0.1 x 0.06mm/Stainless Steel/11cm



The last 5mm on the inside of the tips are highly polished to help prevent tissue from sticking. Ideal for embryos, nerves, or other adhesive tissue.

General Info

Product Number	11412-11
Category	Forceps & Tweezers > Smooth
EAN	4251909502539



Find 11412-11 on finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Style	#5		
Tipdimensions	0.1 x 0.06 mm		
Tips	Standard		
Tipshape	Straight		