

F · S · T

## #3 Forceps - Standard Tips/Straight/Stainless Steel/12cm



Straight forceps with sturdy coarse tips.

## General Info

|                |                             |
|----------------|-----------------------------|
| Product Number | 11445-12                    |
| Category       | Forceps & Tweezers > Smooth |
| EAN            | 4251909502553               |



Find 11445-12 on  
[finescience.com](https://www.finescience.com)

## Technical Info

|                        |                                 |                 |      |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight        | 15 g                            | Net Weight      | 12 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA   |
| Minimum Order Quantity | 1                               |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |      |
| Autoclavesafe          | Yes                             |                 |      |
| Length                 | 12 cm                           |                 |      |
| Tipdimensions          | 0.2 x 0.12 mm                   |                 |      |
| Tips                   | Standard                        |                 |      |
| Tipshape               | Straight                        |                 |      |