



#7 Self-Closing Forceps - Standard Tips/Curved/Stainless Steel/11.5cm



Ideal for lifting and holding EM grids. Consistent pressure at tips. Correctly tensioned crossover shanks come to a microscopically matched point.

General Info

Product Number	11487-11
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	4251909502584



Find 11487-11 on finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	9 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11.5 cm		
Style	#7		
Tipdimensions	0.2 x 0.12 mm		
Tips	Standard		
Tipshape	Curved		