



High precision cutting tweezers are suitable for applications which require cutting small material or performing small cuts. The strong blades are designed to apply moderate cutting force. Suitable for cut small, soft metal strands (i.e. Copper).

General Info

Product Number	11510-12
Category	Forceps & Tweezers > Cutting MedTech Line
EAN	4251909512170



Find 11510-12 on
finescience.com

Technical Info

Shipping Weight	29 g	Net Weight	26 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Carbon Steel (C)		
Autoclavesafe	No		
Bladematerial	Carbon steel		
Cuttingedge	12 mm		
Length	12 cm		
Esdsafe	No		
Tipdimensions	12 x 4.25 mm		
Tips	Cutting		
Tipshape	Angled		