



The high-precision cutting tweezers type #15AGW are suitable for applications requiring the handling of small materials or making fine cuts. The robust blades are designed for applying moderate cutting forces. They are suitable for cutting small, soft metal wires (e.g., copper). The epoxy resin coating gives the tweezers ESD properties, making them ideal for handling electronic components and for use in cleanrooms.

General Info

Product Number	11521-11
Category	Forceps & Tweezers > Cutting MedTech Line
EAN	4251909512194



Find 11521-11 on
finescience.com

Technical Info

Shipping Weight	28 g	Net Weight	25 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Carbon Steel (C), ESD epoxy polyester (NE)		
Autoclavesafe	No		
Bladematerial	Carbon steel		
Cuttingedge	10 mm		
Length	11.5 cm		
Esdsafe	Yes		
Tipdimensions	9.6 x 4 mm		
Tips	Cutting		
Tipshape	Angled		