



#15AGW M Diamond Coated Cutting Tweezers



The diamond-coated cutting forceps type #15AGWM are suitable for applications requiring the cutting of surgical alloys or hard wires in medical technology. The robust blades are designed for the application of moderate cutting forces. The carbon steel ensures hard and stable cutting edges. The blades have a hardness of 61 HRC. The diamond coating (DG4) provides the carbon steel blades with a smooth surface combined with extremely high corrosion resistance as well as good wear and abrasion resistance.

General Info

Product Number	11530-11
Category	Forceps & Tweezers > Cutting MedTech Line
EAN	4251909512224



Find 11530-11 on
finescience.com

Technical Info

Shipping Weight	28 g	Net Weight	25 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Carbon Steel (C), Diamond coating (DG4)		
Autoclavesafe	No		
Bladematerial	Carbon steel		
Cuttingedge	10 mm		
Length	11.5 cm		
Esdsafe	No		
Tipdimensions	9.6 x 4 mm		
Tips	Cutting		
Tipshape	Angled		