



Anatomical Narrow Forceps Titanium - Serrated/Straight/14.5cm



MOST
POPULAR



TITANIUM

Lighter spring action than standard pattern for grasping more delicate tissue.

General Info

Product Number	11602-14
Category	Forceps & Tweezers > Serrated Most Popular > Forceps & Tweezers
EAN	4251909502614



Find 11602-14 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	15 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Titanium		
Antimagnetic	Yes		
Length	14.5 cm		
Tipdimensions	2.25 x 1.45 mm		
Tipwidth	2.25 mm		
Tips	Serrated		
Tipshape	Straight		