



Adson DeBakey Forceps - Atraumatic Serrations/Straight/15cm



Ideal for gentle, atraumatic gripping of vessels, tissue and other delicate structures. Two rows of atraumatic teeth interfit with one opposite row of atraumatic teeth.

General Info

Product Number	11616-15
Category	Forceps & Tweezers > Serrated
EAN	4251909502645



Find 11616-15 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	28 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15 cm		
Tipdimensions	1.75 x 1.5 mm		
Tipwidth	1.5 mm		
Tips	DeBakey serrations		
Tipshape	Straight		