



Adson DeBakey Forceps - Atraumatic Serrations/Straight/12.5cm



Ideal for gentle, atraumatic gripping of vessels, tissue and other delicate structures. Two rows of atraumatic teeth interfit with one opposite row of atraumatic teeth.

General Info

Product Number	11617-12
Category	Forceps & Tweezers > Serrated
EAN	4251909502652



Find 11617-12 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	18 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	12.5 cm		
Tipdimensions	1.75 x 1.5 mm		
Tipwidth	1.5 mm		
Tips	DeBakey serrations		
Tipshape	Straight		