



Dumont #5 Forceps - Standard Tips/Straight/Teflon Coated/11cm



Allgemein

Produktnummer	11626-11
Kategorie	Forceps & Tweezers > Smooth MedTech Line
EAN	7640241754244



Finde 11626-11 auf finescience.com

Technische Informationen

Versandgewicht	15 g	Nettogewicht	12 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Teflon coated Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	11 cm		
Style	#5		
Tipdimensions	0.1 x 0.06 mm		
Tips	Standard		
Tipshape	Straight		