



Dumont #5 Forceps - Standard Tips/Straight/Teflon Coated/11cm



General Info

Product Number	11626-11
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	7640241754244



Find 11626-11 on finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Teflon coated Surgical Stainless Steel (Inox)		
Length	11 cm		
Style	#5		
Tipdimensions	0.1 x 0.06 mm		
Tips	Standard		
Tipshape	Straight		