



Adson-Brown Forceps - Shark Teeth/Straight/12cm



Allgemein

Produktnummer	11627-12
Kategorie	Forceps & Tweezers > Toothed
EAN	4251909502669



Finde 11627-12 auf finescience.com

Technische Informationen

Versandgewicht	25 g	Nettogewicht	21 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	12 cm		
Tipdimensions	2 x 1.4 mm		
Tipwidth	1.4 mm		
Tips	Shark teeth		
Tipshape	Straight		