



Graefe Forceps Titanium - Serrated/Slightly curved/10cm



Allgemein

Produktnummer	11651-10
Kategorie	Forceps & Tweezers > Serrated Most Popular > Forceps & Tweezers
EAN	4251909502683



Finde 11651-10 auf finescience.com

Technische Informationen

Versandgewicht	10 g	Nettogewicht	6 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Titanium		
Antimagnetic	Yes		
Length	10 cm		
Tipdimensions	1.2 x 0.9 mm		
Tipwidth	1.2 mm		
Tips	Serrated		
Tipshape	Slightly curved		