



Graefe Forceps Titanium - Serrated/Curved/10cm



General Info

Product Number	11652-10
Category	Forceps & Tweezers > Serrated Most Popular > Forceps & Tweezers
EAN	4251909502690



Find 11652-10 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	6 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Titanium		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Length	10 cm		
Tipdimensions	1.2 x 0.9 mm		
Tipwidth	1.2 mm		
Tips	Serrated		
Tipshape	Curved		