



Plastic Forceps - Smooth/Sharp Curved/11.5cm



These forceps are acid resistant and able to withstand temperatures up to 464F (240C). Made of acetyl, a strong, stiff plastic characterized by a good fatigue life, low moisture sensitivity, with a high resistance to solvents and chemicals.

General Info

Product Number	11700-01
Category	Forceps & Tweezers > Smooth
EAN	4251909502713



Find 11700-01 on finescience.com

Technical Info

Shipping Weight	70 g	Net Weight	65 g
Package Size	10 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Acetyl		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Length	11.5 cm		
Tipdimensions	0.3 x 0.2 mm		
Tips	Sharp		
Tipshape	Curved		