



Plastic Forceps - Smooth/Blunt Curved/11.5cm



These forceps are acid resistant and able to withstand temperatures up to 464F (240C). Made of acetyl, a strong, stiff plastic characterized by a good fatigue life, low moisture sensitivity, with a high resistance to solvents and chemicals.

General Info

Product Number	11700-04
Category	Forceps & Tweezers > Smooth
EAN	4251909502744



Find 11700-04 on finescience.com

Technical Info

Shipping Weight	70 g	Net Weight	65 g
Package Size	10 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Acetyl		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Length	11.5 cm		
Tipdimensions	1.5 x 0.4 mm		
Tips	Blunt		
Tipshape	Curved		