



An evolutionary series of ergonomic instruments that are lighter than titanium but still have all the benefits of stainless steel. The fully fenestrated spring body is designed with a stability rod to provide balance in your hand.

This product is discontinued and has limited availability.

General Info

Product Number	11800-01
Category	Forceps & Tweezers > Smooth
EAN	4251909509019



Find 11800-01 on
finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	15 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15 cm		
Tipdimensions	0.05 x 0.05 mm		
Tips	Smooth		
Tipshape	Curved		